


## Product Datasheet

### beta IV Tubulin (TUBB4A) antibody (orb1336209)

<b>Description</b>	Mouse monoclonal antibody to beta IV Tubulin (TUBB4A)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Canine, Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF, IHC, WB
<b>Immunogen</b>	Full length human recombinant protein of human TUBB4 (NP_006078) produced in HEK293T cell.
<b>Target</b>	TUBB4A
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>MW</b>	49.4 kDa
<b>Uniprot ID</b>	<a href="#">P04350</a>
<b>NCBI</b>	<a href="#">NM_006087</a>
<b>Dilution Range</b>	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...