


Product Datasheet

DEF8 antibody (orb1336461)

Description	Mouse monoclonal antibody to DEF8
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human DEF8 (NP_997397) produced in HEK293T cell.
Target	DEF8
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	58.5 kDa
Uniprot ID	Q6ZN54
NCBI	NM_207514
Dilution Range	WB 1:500, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...