

## Product Datasheet

### Annexin A3 (ANXA3) antibody (orb1336945)

**Description** Mouse monoclonal antibody to Annexin A3 (ANXA3)

**Species/Host** Mouse

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** FC, IHC, WB

**Immunogen** Full length human recombinant protein of human ANXA3 (NP\_005130) produced in HEK293T cell.

**Target** ANXA3

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Application notes** For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

**Isotype** IgG1

**Clonality** Monoclonal


**MW** 36.2 kDa

**Uniprot ID** [P12429](#)

**NCBI** [NM\\_005139](#)

**Dilution Range** WB 1:500~2000, IHC 1:150, FLOW 1:100

**Expiration Date** 12 months from date of receipt.



Western blot analysis of HEK293T cells w...