

Product Datasheet

Syntrophin alpha 1 (SNTA1) antibody (orb1337003)

Description Mouse monoclonal antibody to Syntrophin alpha 1 (SNTA1)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human SNTA1 (NP_003089) produced in HEK293T cell.

Target SNTA1

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 53.7 kDa

Uniprot ID [Q13424](#)

NCBI [NM_003098](#)

Dilution Range WB 1:2000, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...