


Product Datasheet

PBP (PEBP1) antibody (orb1337386)

Description	Mouse monoclonal antibody to PBP (PEBP1)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PEBP1 (NP_002558) produced in HEK293T cell.
Target	PEBP1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	21.0 kDa
Uniprot ID	P30086
NCBI	NM_002567
Dilution Range	WB 1:2000 IF 1:250
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...