


Product Datasheet

FAHD2A antibody (orb1337982)

Description	Mouse monoclonal antibody to FAHD2A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human FAHD2A (NP_057128) produced in HEK293T cell.
Target	FAHD2A
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	34.6 kDa
Uniprot ID	Q96GK7
NCBI	NM_016044
Dilution Range	WB 1:1000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...