


Product Datasheet

GRASP1 (GRIPAP1) antibody (orb1338077)

| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody to GRASP1 (GRIPAP1) |
| Species/Host | Mouse |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | FC, IF, IHC, WB |
| Immunogen | Full length human recombinant protein of human GRIPAP1 (NP_064522) produced in HEK293T cell. |
| Target | GRIPAP1 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 96.0 kDa |
| Uniprot ID | Q4V328 |
| NCBI | NM_020137 |
| Dilution Range | WB 1:1000, IHC 1:50, IF 1:100, Flow 1:100 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293T cells w...