


Product Datasheet

ERAB (HSD17B10) antibody (orb1338135)

Description	Mouse monoclonal antibody to ERAB (HSD17B10)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human HSD17B10 (NP_004484) produced in HEK293T cell.
Target	HSD17B10
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	26.9 kDa
Uniprot ID	Q99714
NCBI	NM_004493
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...