


Product Datasheet

MAN1 (LEMD3) antibody (orb1338152)

Description	Mouse monoclonal antibody to MAN1 (LEMD3)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IP, WB
Immunogen	Full length human recombinant protein of human LEMD3 (NP_055134) produced in HEK293T cell.
Target	LEMD3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	100.0 kDa
Uniprot ID	Q9Y2U8
NCBI	NM_014319
Dilution Range	WB 1:500, IF 1:100, Flow 1:100, IP 2-4ug/mg
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...