
Product Datasheet


Protein Kinase A regulatory subunit I alpha (PRKAR1A) antibody (orb1338292)

Description

Mouse monoclonal antibody to Protein Kinase


Species/Host

Mouse


Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IF, IHC, WB

Western blot analysis of HEK293T cells w...

Immunogen

Full-length protein expressed in 293T cell transfected with human PRKAR1A expression vector

Target

PRKAR1A

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype

IgG1

Clonality

Monoclonal

MW

43.0 kDa

Uniprot ID
[P10644](#)
NCBI
[NM_212472](#)
Dilution Range

WB 1:2000, IHC 1:50, IF 1:100

Expiration Date

12 months from date of receipt.