


Product Datasheet

TrkC (NTRK3) antibody (orb1338392)

Description	Mouse monoclonal antibody to TrkC (NTRK3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human NTRK3 expression vector
Target	NTRK3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	94.4 kDa
Uniprot ID	Q16288
NCBI	NM_001007156
Dilution Range	WB 1:1000 IHC 1:50 IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...