


Product Datasheet

NEK6 antibody (orb1338401)

Description	Mouse monoclonal antibody to NEK6
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human NEK6 (NP_055212) produced in HEK293T cell.
Target	NEK6
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	35.5 kDa
Uniprot ID	Q9HC98
NCBI	NM_014397
Dilution Range	WB 1:1000, IHC 1:50, IF 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...