


Product Datasheet

Grp75 (HSPA9) antibody (orb1338413)

Description	Mouse monoclonal antibody to Grp75 (HSPA9)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human HSPA9 (NP_004125) produced in HEK293T cell.
Target	HSPA9
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	73.5 kDa
Uniprot ID	P38646
NCBI	NM_004134
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:50, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...