


Product Datasheet

NEK6 antibody (orb1338418)

Description	Mouse monoclonal antibody to NEK6
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human NEK6 expression vector
Target	NEK6
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	35.7 kDa
Uniprot ID	Q9HC98
NCBI	NM_014397
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, Flow 1:500, IP 2-4ug/mg
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...