


Product Datasheet

Kindlin 2 (FERMT2) antibody (orb1338424)

Description	Mouse monoclonal antibody to Kindlin 2 (FERMT2)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human FERMT2 expression vector
Target	FERMT2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	77.9 kDa
Uniprot ID	Q96AC1
NCBI	NM_006832
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:50~100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...