


Product Datasheet

Rad9 (RAD9A) antibody (orb1338490)

Description	Mouse monoclonal antibody to Rad9 (RAD9A)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human RAD9A expression vector
Target	RAD9A
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	42.5 kDa
Uniprot ID	Q99638
NCBI	NM_004584
Dilution Range	WB 1:2000 IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...