


Product Datasheet

TIF1 gamma (TRIM33) antibody (orb1338634)

Description	Mouse monoclonal antibody to TIF1 gamma (TRIM33)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 887-1046 of human Trim33 (NP_056990) produced in E.coli.
Target	TRIM33
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	122.5 kDa
Uniprot ID	Q9UPN9
NCBI	NM_015906
Dilution Range	WB 1:500, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...