


Product Datasheet

Sonic Hedgehog (SHH) antibody (orb1338678)

Description	Mouse monoclonal antibody to Sonic Hedgehog (SHH)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 24-197 of human SHH
Target	SHH
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	27.6 kDa (Actual 22 kDa/55 kDa)
Uniprot ID	Q15465
NCBI	NM_000193
Dilution Range	WB: 1:2000, IF (1:100)
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...