


Product Datasheet

Hex (HHEX) antibody (orb1338691)

Description	Mouse monoclonal antibody to Hex (HHEX)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fragment expressed in E.coli corresponding to amino acids 100-270 of human Hex
Target	HHEX
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	29.8 kDa (experimental 44 kDa)
Uniprot ID	Q03014
NCBI	NM_002729
Dilution Range	WB: 1:50 - 1:100, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...