


Product Datasheet

TAG-72 antibody (orb1338782)

Description	Mouse monoclonal antibody to TAG-72
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Membrane-enriched fraction of a human breast carcinoma liver metastasis
Target	TAG-72
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...