
Product Datasheet

Envelope Glycoprotein Yellow Fever Virus antibody (orb13392)

Description

Envelope Glycoprotein Yellow Fever

Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from Yellow fever virus envelope glycoprotein E (601-700/3411 aa)
Target	Yellow fever virus envelope glycoprotein E
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54/375 kDa
Uniprot ID	P03314
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.