


Product Datasheet

RASGRP 4 (RASGRP4) Human Over-expression Lysate (orb1339776)

| | |
|------------------------|--|
| Description | Transient overexpression lysate of RAS guanyl releasing protein 4 (RASGRP4), transcript variant d |
| Target | RASGRP4 |
| Storage | The lysate is shipped with dry ice. The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C. |
| Tag | C-Myc/DDK |
| Note | For research use only |
| Source | Human |
| MW | 66.8 kDa |
| Uniprot ID | COLTP7 , Q8TDF6 |
| Expiration Date | 6 months from date of receipt. |



Western blot analysis of Left: Cell lysa...