


Product Datasheet

TCF7L2 Human Over-expression Lysate (orb1339778)

Description	Transient overexpression lysate of transcription factor 7-like 2 (T-cell specific, HMG-box) (TCF7L2)...
Target	TCF7L2
Storage	The lysate is shipped with dry ice. The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	65.9 kDa
Uniprot ID	Q9NQB0
Expiration Date	6 months from date of receipt.



Western blot analysis of
Left: Cell lysa...