

## Product Datasheet

### Endothelin 1 antibody (orb13403)

**Description**

Endothelin 1 antibody

**Species/Host**

Rabbit

**Reactivity**

Bovine, Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat, Sheep

**Conjugation**

Unconjugated

**Tested**

ELISA, ICC, IF, WB

**Applications**

**Immunogen**

KLH conjugated synthetic peptide derived from human Endothelin 1. Please contact us for the exact immunogen sequence. The peptide is available as orb374764.

**Target**

Endothelin 1

**Form/Appearance**

10 mM PBS, 0.02% sodium azide

**Concentration**

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

**Storage**

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Polyclonal antibodies are purified by peptide affinity chromatography

**MW**

24 kDa

**Uniprot ID**

**P05305, P22388, P22387**

**NCBI**

**001168319, 001161791**

**Entrez**


**1906**

**Dilution Range**


IF/ICC: 1:200, WB: 1:100-1000

**Expiration Date**


12 months from date of receipt.



Western blot analysis of Hela cells (lan...



Immunofluorescence image of human hela c...



IF analysis of human hela cell tissue us...