


## Product Datasheet

# PIP5K1 alpha (PIP5K1A) Human Over-expression Lysate (orb1340787)

<b>Description</b>	Transient overexpression lysate of phosphatidylinositol-4-phosphate 5-kinase, type I, alpha (PIP5K1A...)
<b>Target</b>	PIP5K1A
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q99755</a>
<b>Expiration Date</b>	6 months from date of receipt.



Western blot analysis of Left: Cell lysa...