


Product Datasheet

OR4D11 Human Over-expression Lysate (orb1342247)

Description	Transient overexpression lysate of olfactory receptor, family 4, subfamily D, member 11 (OR4D11)
Target	OR4D11
Storage	The lysate is shipped with dry ice. The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	34.8 kDa
Uniprot ID	Q8NGI4 , A0A126GVQ9
Expiration Date	6 months from date of receipt.



Western blot analysis of Left: Cell lysa...