


Product Datasheet

APPBP1 (NAE1) Human Over-expression Lysate (orb1342937)

Description	Transient overexpression lysate of NEDD8 activating enzyme E1 subunit 1 (NAE1), transcript variant 3
Target	APPBP1
Storage	Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	60.2 kDa
Uniprot ID	Q13564
NCBI	NP_001018170, NM_001018160
Expiration Date	6 months from date of receipt.



Western blot analysis of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using APPBP1 (NAE1) Human Over-expression Lysate

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)