

Product Datasheet

GPR109A antibody (orb13445)

Description

GPR109A antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from GPR109A. Please contact us for the exact immunogen sequence. The peptide is available as orb393919.

Target

GPR109A

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

42 kDa

Uniprot ID

[Q8TDS4](#)

Entrez


[338442](#)

Dilution Range


IHC-P: 1:100-500 (based on 0.5 mg/ml), WB: 1:100-500

Expiration Date

12 months from date of receipt.




Immunohistochemical staining of rat lymph...



WB image of paraffin embedded rat lymph

...



IHC-P image of rat lymph node tissue using...