


Product Datasheet

RAC2 Human Over-expression Lysate (orb1345226)

Description	Transient overexpression lysate of ras-related C3 botulinum toxin substrate 2 (rho family, small GTP...
Target	RAC2
Storage	The lysate is shipped with dry ice. The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	2.6 kDa
Uniprot ID	V9H0H7 , P15153 , A0A024R1P2
Expiration Date	6 months from date of receipt.



Western blot analysis of
Left: Cell lysa...