


## Product Datasheet

# SPT3 (SUPT3H) Human Over-expression Lysate (orb1345590)

<b>Description</b>	Transient overexpression lysate of suppressor of Ty 3 homolog ( <i>S. cerevisiae</i> ) (SUPT3H), transcript v...
<b>Target</b>	SUPT3H
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>A0A024RD67, O75486</b>
<b>Expiration Date</b>	6 months from date of receipt.



Western blot analysis of Left: Cell lysa...