


## Product Datasheet

# PPFIA3 Human Over-expression Lysate (orb1345636)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Transient overexpression lysate of protein tyrosine phosphatase, receptor type, f polypeptide (PTPRF... |
| <b>Target</b>          | PPFIA3  |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>075145</b>   |
| <b>Expiration Date</b> | 6 months from date of receipt.  |



Western  
blot  
analysis of  
Left: Cell  
lysa...