


## Product Datasheet

# HPS2 (AP3B1) Human Over-expression Lysate (orb1345639)

<b>Description</b>	Transient overexpression lysate of adaptor-related protein complex 3, beta 1 subunit (AP3B1)
<b>Target</b>	AP3B1
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">A0A0S2Z5J4</a> , <a href="#">O00203</a>
<b>Expiration Date</b>	6 months from date of receipt.



Western blot analysis of  
Left: Cell lysa...