


## Product Datasheet

# Cyclophilin F (PPIF) Human Over-expression Lysate (orb1346668)

<b>Description</b>	Transient overexpression lysate of peptidylprolyl isomerase F (PPIF), nuclear gene encoding mitochon...
<b>Target</b>	PPIF
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">A0A024QZS4</a> , <a href="#">P30405</a>
<b>Expiration Date</b>	6 months from date of receipt.



Western blot analysis of Left: Cell lysa...