

Product Datasheet

POLD2 antibody (orb134912)

Description	Rabbit polyclonal to POLD2. The function of the small subunit is not yet clear.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	POLD2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51289
Uniprot ID	P49005
NCBI	NP_001120690.1 , NP_001243808.1 , NP_006221.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

