

Product Datasheet

BRTC1/2 antibody (orb134924)

Description	Rabbit polyclonal to BRTC1/2. Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	BTRC
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	68867
Uniprot ID	Q9Y297
NCBI	NP_003930.1 , NP_378663.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

