



Product Datasheet

SUMO2 antibody (orb134982)

Description	Rabbit polyclonal to SUMO2. Ubiquitin-like protein that can be covalently attached to proteins as a ...
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant Protein
Target	SUMO2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	10871
Uniprot ID	P61956
NCBI	NP_001005849.1 , NP_008868.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...