

Product Datasheet

NEK3 Human Over-expression Lysate (orb1353508)

Description	Transient overexpression lysate of NIMA (never in mitosis gene a)-related kinase 3 (NEK3), transcrip...
Target	NEK3
Storage	The lysate is shipped with dry ice. The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	57.7 kDa
Uniprot ID	Q6ZN64 , P51956
Expiration Date	6 months from date of receipt.