

## Product Datasheet

# UBE2D2 Human Over-expression Lysate (orb1354639)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Transient overexpression lysate of ubiquitin-conjugating enzyme E2D 2 (UBC4/5 homolog, yeast) (UBE2D2), transcript variant 2   |
| <b>Target</b>          | UBE2D2   |
| <b>Storage</b>         | Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C. |
| <b>Tag</b>             | C-Myc/DDK  |
| <b>Note</b>            | For research use only  |
| <b>Source</b>          | Human  |
| <b>MW</b>              | 13.5 kDa   |
| <b>Uniprot ID</b>      | <b>P62837</b>  |
| <b>NCBI</b>            | <b>NM_181838, NP_862821</b>  |
| <b>Expiration Date</b> | 6 months from date of receipt.   |

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)