

Product Datasheet

SFMBT1 Human Over-expression Lysate (orb1356792)

Description	Transient overexpression lysate of Scm-like with four mbt domains 1 (SFMBT1), transcript variant 3
Target	SFMBT1
Storage	The lysate is shipped with dry ice. The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	98 kDa
Uniprot ID	A0A024R338 , Q9UHJ3
Expiration Date	6 months from date of receipt.