



Product Datasheet

LAMP1 antibody (orb135825)


Description	Rabbit polyclonal antibody to LAMP1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human LAMP1 (NP_005552.3).
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	90-120kDa
Uniprot ID	P11279
Entrez	3916
Dilution Range	WB: 1:500-,2000 IF: 1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of NIH/3T3 C...



Western blot analysis of extracts of mou...