


Product Datasheet

TCF7L2 Human Over-expression Lysate (orb1359191)

Description	TCF7L2 HEK293T cell transient overexpression lysate (as WB positive control)
Target	TCF7L2
Storage	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	45.8 kDa
Uniprot ID	C6ZRK2, A0A0A0MTL7
Expiration Date	6 months from date of receipt.



Western blot validation of Left: Cell ly...