


## Product Datasheet

# CT45A5 Human Over-expression Lysate (orb1359631)

<b>Description</b>	CT45A5 HEK293T cell transient overexpression lysate (as WB positive control)
<b>Target</b>	CT45A5
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q6NSH3</a> , <a href="#">P0DMV0</a> , <a href="#">P0DMU7</a> , <a href="#">P0DMU8</a>
<b>Expiration Date</b>	6 months from date of receipt.



Western blot validation of Left: Cell ly...