


## Product Datasheet

# WDR21A (DCAF4) Human Over-expression Lysate (orb1359733)

<b>Description</b>	DCAF4 HEK293T cell transient overexpression lysate (as WB positive control)
<b>Target</b>	DCAF4
<b>Storage</b>	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
<b>Tag</b>	C-Myc/DDK
<b>Note</b>	For research use only
<b>Source</b>	Human
<b>MW</b>	48.6 kDa
<b>Uniprot ID</b>	<a href="#">Q8WV16</a>
<b>Expiration Date</b>	6 months from date of receipt.



Western blot validation of Left: Cell ly...