

Product Datasheet

SAE1 Human Over-expression Lysate (orb1360347)

Description	SAE1 HEK293T cell transient overexpression lysate (as WB positive control)
Target	SAE1
Storage	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	38.4 kDa
Uniprot ID	Q9UBE0
Expiration Date	6 months from date of receipt.