



## Product Datasheet

### RBP4 antibody (orb136234)


<b>Description</b>	Rabbit polyclonal antibody to RBP4
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant fusion protein of human RBP4 (NP_006735.2).
<b>Preservatives</b>	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>MW</b>	23kDa
<b>Uniprot ID</b>	<b>P02753</b>
<b>Entrez</b>	<b>5950</b>
<b>Dilution Range</b>	WB: 1:500-2000, IHC-P: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.



WB analysis of RBP4 (left) and antibody...



Immunohistochemical analysis of formalin...



Western blot analysis of extracts of var...