


Product Datasheet

OR2T5 Human Over-expression Lysate (orb1362591)

| | |
|------------------------|---|
| Description | OR2T5 HEK293T cell transient overexpression lysate (as WB positive control) |
| Target | OR2T5 |
| Storage | The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C. |
| Tag | C-Myc/DDK |
| Note | For research use only |
| Source | Human |
| MW | 35.4 kDa |
| Uniprot ID | Q6IEZ7 |
| Expiration Date | 6 months from date of receipt. |



Western blot validation of Left: Cell ly...