



Product Datasheet

SSTR2 antibody (orb136360)


Description	Rabbit polyclonal antibody to SSTR2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human Somatostatin Receptor 2 (SSTR2)(NP_001041.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	87KDa
Uniprot ID	P30874
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human li...



Western blot analysis of extracts of mou...



Immunohistochemical staining of human st...