


## Product Datasheet

# IRAG1 Human Over-expression Lysate (orb1364485)

|                        |   |
|------------------------|---|
| <b>Description</b>     | MRV11 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | IRAG1   |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <a href="#">Q9Y6F6</a>  |
| <b>Expiration Date</b> | 6 months from date of receipt.  |



Western blot validation of Left: Cell ly...