

Product Datasheet

NEK3 Human Over-expression Lysate (orb1365374)

Description	NEK3 HEK293T cell transient overexpression lysate (as WB positive control)
Target	NEK3
Storage	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
Tag	C-Myc/DDK
Note	For research use only
Source	Human
MW	57.5 kDa
Uniprot ID	P51956
Expiration Date	6 months from date of receipt.