


## Product Datasheet

# RAC2 Human Over-expression Lysate (orb1365566)

<b>Description</b>	RAC2 HEK293T cell transient overexpression lysate (as WB positive control)
<b>Target</b>	RAC2
<b>Storage</b>	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
<b>Tag</b>	C-Myc/DDK
<b>Note</b>	For research use only
<b>Source</b>	Human
<b>MW</b>	2.6 kDa
<b>Uniprot ID</b>	<a href="#">V9H0H7</a> , <a href="#">P15153</a> , <a href="#">A0A024R1P2</a>
<b>Expiration Date</b>	6 months from date of receipt.



Western blot validation of Left: Cell ly...